

COMPRESSION STOCKING **AD LIPOTROMB**

Worried about thrombosis after
liposuction of thighs?

No need for that
with **AD LIPOTROMB**
compression
stockings!



COMPRESSION
18 – 21 mmHg

For ages, LIPOELASTIC has been taking care of **prevention during plastic surgery interventions.**

LIPOTROMB thigh stockings are already a well established prevention of thrombosis for patients who are immobile. However their use was always problematic for patients who had undergone liposuction in the area of the knees and thighs. Which is why we've developed AD LIPOTROMB compression calf stockings.

- for prevention of thrombosis with the correct medical compression, tightest at the ankle, with compression gradually decreasing up the leg
- enhances the safety of the surgical intervention
- the compression creates pressure on the skin surface and deep veins in the legs, thus regulating blood flow towards the heart and keeping circulation from stagnating
- reduces the occurrence of leg swelling
- number range weaved into fabric and letter system of size labeling enables the correct identification of the product
- breathable and pleasant to skin
- made from material resistant to oils and ointments

Color: white



full tip with an inspection opening under toes



finely ribbed upper edge



size weaved into fabric and reinforced material in the heel area

Size	S	M	L	XL	XXL	XXXL	4XL
d	31 – 38	33 – 41	35 – 44	37 – 47	39 – 49	41 – 51	52 – 62
b	21 – 23	23 – 25	25 – 27	27 – 29	29 – 31	30 – 33	34 – 38

d – measure the calf circumference in the place where the edge is supposed to be

b – measure the ankle circumference in the narrowest place



LIPOELASTIC a.s., Radlická 608/2, 150 00 Praha 5, CZECH REPUBLIC
E-mail: info@lipoelastic.com | www.lipoelastic.com | MADE IN EU